# SinCrest

technique

CE  
0123

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### Intended use

SinCrest is a kit for the maxillary sinus lift with crestal approach. The SinCrest system consists of a kit including handpiece drills, depth stops at different levels and a manual osteotome. The manual osteotome was designed to obtain the controlled fracture of the bone floor near the maxillary sinus membrane without damaging the membrane itself. Thus dental implants can be inserted exceeding the initial bone length.

- **SinCrest Kit is not sterile when supplied. It must be washed and sterilised every time before using it.**
- **The device must be used exclusively by skilled medical staff only.**
- **The surgeon must assess whether the patient is suitable for osteotomy surgery and choose the correct surgical sequence.**
- **Never use SinCrest without making sure that all its components have been assembled correctly (see figure n. 16)**
- **Before using SinCrest, press the probe onto a sterile field and check the correct recovery by means of the spring. (see figure n. 1)**
- **Make sure that a dental floss is inserted onto the RX Pins to avoid the risk of accidental ingestion by the patient. o An incorrect speed (exceeding 850 RPM) of the hand-piece may overheat the bone tissue leading to bone necrosis.**
- **Attention, make sure that the Drill Gripper device is not reassembled during the sterilization stage, but only before using it (see figure n. 15).**
- **META does not shoulder any responsibility for the incorrect use of the product.**


Fig. 1

### Contraindications

- Do not use the SinCrest technique in case of residual bone shorter than 5mm and taller than 11mm.
- We recommend not to use the SinCrest device in case of reduced stability during the screwing stage in the bone hole drilled
- Do not use the SinCrest device if the vertical distance from the maxillary sinus floor to the terminal occlusion of the tooth adjacent to the site to be treated exceeds 23mm (see figure n. 2).
- The Ref 4344 Probe Drill Ø 1.2 is not designed to cut the bone laterally and must be used only to drill in depth. An incorrect use may damage the drill.


Fig. 2

- Do not use the SinCrest device without having previously drilled a hole with the relevant flat SinCrest Drills Ø 3.0 or SinCrest Drill Ø 4.0 (see figure n. 3)


Fig. 3

- Do not use tools such as pliers or wrenches to screw the SinCrest device into the hole in the bone (see figure n. 4)

Fig. 4

## SINCREST KIT USAGE PROCEDURE

### Precautions

Correct X-rays, direct palpation and inspections of the operating site are necessary to identify the anatomy of the available bone. The most suitable surgical sequence is to be assessed according to the available bone and the size of the dental implant. During the use of the drills to prepare the bone hole, ensure an abundant external lavage with sterile saline solution; a good lavage may also be ensured by means of a sterile syringe.

## Surgical procedure

After having proceeded to the full depth incision and elevation of the soft tissues, the site may be prepared to use the selected SinCrest device. To handle the drills, we suggest to use the accessory called "Drill Gripper", which allows preventing the direct contact with the operator's hands. From the black side the Ref 4338-SinCrest drill Ø 3.2, Ref 4381-SinCrest drill Ø 4.0 and Ref 4340-SinCrest drill HB Ø 4.2 may be withdrawn, whereas the brown side contains all the other drills of the SinCrest Kit (see figure n. 5)


Fig. 5

The laser marks reported on the Probe drill Ø 1.2, the SinCrest drills Ø 3.0, 3.2, 4.0 and 4.2 and the RX Pin pins indicate the operating depth at 4-5-6-7-8-9-10 and 11mm. (see figure n. 6)


Fig. 6

**Ref 4343-Locator drill:** use the drill to identify the correct position of the implant and remove the cortical bone tissue at a 3.5mm depth as shown by the laser mark.


**Ref 4344-Probe drill Ø 1.2:** to be used after having positioned the depth stop selected according to the bone availability (see paragraph on Stop positioning). Since the osteotomy achieved defines the final position and angle of the implant, it must be carried out with a view to the full prosthetic plan. Use the drill with back and forth movements until the desired depth is

achieved. Only the head of the device is sharp. No lateral cutting strength can be implemented to modify the osteotomy angle: it must be used only to drill in depth.

**Ref 4354-RX pin Ø 1.2:** use the pins to check the correct angle and depth achieved through an x-ray and/or visual examination.

**Ref 4379-Sin Probe Ø 1.2:** in addition to the tactual check, this device allows a metric assessment of the operating depth.

**Ref 4380-Guide drill Ø 3.0:** use the drill to create a coronal recess in the cortical layer to achieve a 2 mm depth as defined by the stop on the drill. This drill allows the correct centring of the SinCrest Drill (see figure n. 7).

**Ref 4337-SinCrest drill Ø 3.0:** it may be used after having positioned the depth stop selected according to the available bone (see paragraph on Stop positioning). To obtain a non-flared hole which may jeopardize the stabilization of the SinCrest device, operate the motor only after having placed the SinCrest drill in the guide hole (see figure n. 7). Prepare the site until the pre-determined depth is achieved. To obtain a precise hole diameter, use the drill only once until the stop stroke is achieved. Further uses may lead to a poorer stabilization of the SinCrest device, in particular if little vertical bone is available. We recommend to maintain an approximately 2.0 mm distance from the maxillary sinus floor membrane (see figure n. 3).

**Ref 4356-RX pin Ø 3.0:** it confirms the correct angle


and depth achieved during the previous osteotomy stages.

#### Notice:


During this stage the surgeon must decide whether the surgical sequence should be implemented to use the SinCrest Ø 4.0 so as to place a dental implant with greater diameter. If necessary, refer to the Ref 4381 SinCrest Ø 4.0 operating instructions.

#### Ref 4333-SinCrest Ø 3.0 - brown colour


**Step 0** Gradually screw the device inside the dedicated guide hole. The device has achieved the depth obtained by means of the drills when the white mark emerges minimally during the progress of the probe as compared to the grip (see figure n. 8).

**Step 1** Turn the handle anti-clockwise by 1/2 revolution and at the same time implement an axial pressure, start the osteotomy with "back and forth" turning movements by 1/2 turn until the distance (0.5mm) between the tooth and the stroke is zero (see figure n. 9). The SinCrest device has progressed by 0.5 mm in the apical direction.


By manually implementing pressure on the probe (see figure n. 10), constantly check the residual strength of the maxillary sinus floor.

The progress in the apical direction of the probe shows that the bone operculum has been fractured; unscrew the SinCrest device and insert bio-material or directly the dental implant. In case of residual bone floor strength detected by the probe, proceed to Step 2.

**Step 2** Screw by 1/2 turn (each 1/2 turn equals a 0.5 mm progress) after the tooth is engaged in the stroke (see figure n. 11) and repeat the steps 1 and 2 until the fracture of the maxillary sinus floor reached.

**Ref 4338-SinCrest drill HB Ø 3.2:** this drill must be used in case of hard bone encountered during the screwing stage of the SinCrest Ø 3.0 device. This drill does not require positioning a stop since it must simply widen the hole without working in depth. The osteotomy may be resumed by means of the SinCrest Ø 3.0 device. If bone is available, the surgical procedure may be continued as illustrated in the following steps so as to

position a dental implant with a wide diameter.

**Ref 4381-SinCrest drill Ø 4.0:** the use of this drill must always be preceded by the SinCrest drill Ø 3.0 and SinCrest drill HB Ø 3.2. This drill does not require a stop since it must simply widen the hole without working in depth. To obtain a precise hole diameter, use the drill only once until the stop stroke is achieved. Further uses may lead to a poorer stabilization of the SinCrest device, in particular if little vertical bone is available.

**Ref 4357-RX pin Ø 4.0:** confirms the correct angle and depth achieved during osteotomy.

**Ref 4334-SinCrest Ø 4.0 - black colour.**

**See the Ref 4333-SinCrest Ø 3.0**


**Ref 4340-SinCrest drill HB Ø 4.2:** this drill must be used in case of hard bone encountered during the screwing stage of the SinCrest Ø 4.0 device. This drill does not require positioning a stop since it must simply widen the hole without working in depth. Prepare the site until the desired depth is achieved. The osteotomy may be resumed by means of the SinCrest Ø 4.0 device.

## USAGE PROCEDURE SINCREST MINI-KIT

**Ref 4465 SinCrest Mini-Kit Ø 3.0.**

**Ref 4466 SinCrest Mini-Kit Ø 4.0.**

To ensure correct centring, the use of the Ref 4337 SinCrest Drill Ø 3.0 must with flat end and sharp tip must be preceded by a preparation by means of an Ø 3.0mm(see figure n. 12).

Proceed according to the surgical sequence starting from Ref 4337-SinCrest drill HB Ø 3.0

### Stop positioning

8 stops are available for different depths (4-5-6-7-8-9-10 and 11m) with laser marks for easier identification. The stops are organised in sequence in the dedicated housings in the tray. Their positioning on the Probe Drill Ø 1.2 and SinCrest drill Ø 3.0 takes place directly by simply inserting the drill (see figure n. 13) mounted on the contra-angle hand-piece, directly in the desired depth stop.


Fig. 13

The stop already used may be temporarily positioned in the stop housing, thus avoiding the direct contact with the operator's hands. (see figure n. 14) To avoid the accidental fall of the stops, slide the cover into the lock position.

Fig. 14

### Washing and disinfecting

SinCrest Kit is not sterile when supplied. It must be washed and sterilised each time before it is used. Disinfect all the SinCrest Kit devices by means of the dedicated germicide detergents immediately after the surgical procedure to prevent cross contamination. Disassemble the Drill Gripper device to wash its internal parts (see figure n. 15) and wash the inner cavity by with the help of a brush.


Fig. 15

To disassemble the SinCrest Ø 3.0 and Ø 4.0, use the Screwdriver accessory contained in the tray which allows unscrewing the container (see figure n. 16). Remove the PROBE and SPRING internal components and wash the inner cavity by with the help of a brush.


Fig. 16

To wash, use suitable detergents. Do not use nitrogen peroxide ( $H_2O_2$ ) nor chlorine-based detergents. Wash the devices accurately with flowing water to remove the residual tissue with the help of a soft brush and check that they are not damaged. Disinfect and wash the Tray each time immediately after it has been used. Make sure that no residues remain in the drill and tool housings; in certain cases, disassembling the tools is recommended for better cleaning. When the Tray is washed in dedicated tool-washers, place the tray diagonally for a better elimination of the condensate steam.

## **Sterilization and storage (see Device sterilisation and disinfection list)**

Before proceeding with the sterilisation stage, correctly re-assemble the components and make sure that the spring is repositioned as shown in figure 16 (smaller diameter facing the internal part of the device) Make sure that the Drill Gripper device is not re-assembled during the sterilisation, but only before using it. Use the dedicated SinCrest Tray to arrange the surgical devices. Seal the tray or individual components individually in a suitable sterile bag. Sterilisation should take place inside a steam steriliser at 121°C during 30 min. always using validated sterilisation cycles according to the laws in force. The maintenance of sterilisation is ensured by the a whole bag: always check that it is not damaged before using it.

**You can avoid mistakes and ensure the correct operation of our product only if you read the present document.**

**The present operating instructions and additional information on the SinCrest technique are available on our web site: [www.metahosp.com](http://www.metahosp.com)**

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## List of sterilization and disinfection devices

SinCrest Kit Ø 3.0/4.0 Ref 4418	SinCrest Kit Ø 3.0 Ref 4417	Product	Ref	Description	Reusable	Steam sterilisation	Disassembling recommended	Material
●	●		4343	LOCATOR DRILL	Yes <sup>1,2</sup>	Yes		Stainless Steel
●	●		4344	PROBE DRILL Ø 1.2	Yes <sup>1,2,3</sup>	Yes		Stainless Steel
●	●		4380	GUIDE DRILL Ø 3.0	Yes <sup>1,2</sup>	Yes		Stainless Steel
●	●		4337	SINCREST DRILL Ø 3.0	Yes <sup>1,2,3</sup>	Yes		Stainless Steel
●	●		4381	SINCREST DRILL Ø 4.0	Yes <sup>1,2</sup>	Yes		Stainless Steel
●	●		4338	SINCREST DRILL HB Ø 3.2	Yes <sup>1,2</sup>	Yes		Stainless Steel
●	●		4340	SINCREST DRILL HB Ø 4.2	Yes <sup>1,2</sup>	Yes		Stainless Steel
●	●		4354	RX PIN Ø 1.2	Yes	Yes		Ti6Al4V
●	●		4356	RX PIN Ø 3	Yes	Yes		Ti6Al4V
●	●		4357	RX PIN Ø 4	Yes	Yes		Ti6Al4V
●	●		4379	SIN PROBE Ø1.2	Yes	Yes		Ti6Al4V
●	●		4333	SINCREST Ø 3.0	Yes <sup>4</sup>	Yes	Yes	Ti6Al4V/ Stainless Steel/ Polymer
●	●		4334	SINCREST Ø 4.0	Yes <sup>4</sup>	Yes	Yes	Ti6Al4V/ Stainless Steel/ Polymer
●	●		4339	DRILL GRIPPER	Yes	Yes	Yes	Polymer
●	●		4341	SCREWDRIVER	Yes	Yes		Polymer
●	●		4406 4407 4408 4409 4410 4411 4412 4413	DRILL STOP 4-5-6-7-8-9-10-11mm	Yes	Yes		Titanium
●	●		4439	SINCREST STOP TRAY	Yes	Yes	Yes	Polymer/ Stainless Steel
●	●		4342	SINCREST TRAY	Yes	Yes	Yes	Polymer/ Stainless Steel

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